# MAKRITE SEKURA-N95


## Package :

1. Size : L415×W300×H385mm
2. Volume : 0.0479cbm
3. Quantity : 40PCS×12Box
4. GW : 4.40KG

# Products Introduction

## Packing bag

1. Quantity: 40pcs/bag
2. Weight: 0.25Kg
3. Size: 150×95mm

## The packaging box

1. Quantity: 40pcs/Box
2. Weight: 0.3Kg
3. Size: 142×132×182mm

## The packaging carton

1. Quantity: 12box/ctn
2. Gross Weight: 4.40Kg
3. Size: 415×300×385mm
4. Volume: 0.0479cbm

Product six-sided pictures

**Makrite**<sup>®</sup>  
Life Holds Unlimited Possibilities

*NIOSH  
Approved  
Product*

# SEKURA-N95

## Disposable Respirator N95

95% Efficient Protection Against  
Non-Oil Particulates

**LATEX**

**40**  
Particulate  
Respirators

## RESPIRATOR FITTING INSTRUCTION

1

### STEP 1

Hold the respirator in hand with the nosepiece at your fingertips, allowing the headbands to hang freely below your hand.

2

### STEP 2

Press the respirator firmly against your face with the nosepiece on the bridge of your nose.

3

### STEP 3

Stretch and position the top band high on the back of the head. Stretch the bottom band over the head and position below your ears.

4

### STEP 4

Using both hands, mold nosepiece to the shape of your nose.

5

### STEP 5

To test fit: a) Cup both hands over the respirator being careful not to disturb position, and b) inhale vigorously. If air leaks around the edges, reposition the straps or adjust strap tension for better fit.

Please carefully follow these fitting instructions during each use to achieve proper fit.

MAKRITE INDUSTRIES INC.

11F-5, No.79, Sec. 1, Xintai 5th Rd., XiZhi Dist., New Taipei City 221, Taiwan

Tel: +886-2-2698-2419

Made In China

**Makrite**<sup>®</sup>  
Life Holds Unlimited Possibilities

*NIOSH  
Approved  
Product*

# SEKURA-N95

## Disposable Respirator N95

95% Efficient Protection Against  
Non-Oil Particulates

**40**  
Particulate  
Respirators

MAKRITE INDUSTRIES INC.  
11F-5, No.79, Sec. 1,  
Xintai 5th Rd., XIZhi Dist.,  
New Taipei City 221, Taiwan  
Tel: +886-2-2698-2419

This respirator is approved only in the following configuration:

TC <sup>1</sup>	Protection <sup>1</sup>	Respirator	Cautions and Limitations <sup>2</sup>
		SEKURA-N95	
84A-7638	N95	x	ABCJMNOP

1. Protection

N95-Particulate Filter (95% filter efficiency level) effective against particulate aerosols free of oil; time use restrictions may apply

2. Cautions and Limitations

- A - Not for use in atmospheres containing less than 19.5% oxygen.
- B - Not for use in atmospheres immediately dangerous to life or health.
- C - Do Not exceed maximum use concentrations established by regulatory standards
- J - Failure to properly use and maintain this product could result in injury or death.
- M - All approved respirators shall be selected, fitted, used, and maintained in accordance with MSHA, OSHA and other applicable regulations
- N - Never substitute, modify, add, or omit parts. Use only exact replacement parts in the configuration as specified by the manufacturer.
- O - Refer to users instructions, and/or maintenance manuals for information on use and maintenance of these respirators
- P - NIOSH does not evaluate respirators for use as surgical masks.

### SERVICE TIME LIMITATION

1. Do not use the same respirator for more than 8 hours, continuous or intermittent, in dirty workplaces that could result in high filter loading.
2. Leave contaminated area and remove respirator if breathing becomes difficult or proper fit cannot be obtained or the respirator is damaged, soiled, or distorted. Do not use the same respirator again once it has been used in contaminated area.
3. The maximum service time of this respirator is five (5) consecutive calendar days (including days of non-use), beginning from the first day of use.

### LIMITED WARRANTY IMPORTANT NOTICE TO PURCHASER:

This limited warranty is made in lieu of the warranties of merchantability, fitness for particular purposes and all other warranties, expressed or implied. There are no other warranties which extend beyond the description on the face hereof. The physical standards and specifications of Makrite will be met by products sold.

Exclusive Remedies: damages for the breach of this limited warranty are limited to the replacement of such quantity of Makrite, products proved to be defectively manufactured. Except as provided above, Makrite shall not be liable or responsible for any loss, damage, or liability, direct, indirect, incidental, special, or consequential, arising out of sale, use, or misuse, or the inability to use products by the user.

Made in China

### WARNINGS:

1. Wearer must read fitting instructions and use limitations before use.
2. This respirator is not guaranteed to protect against the risk of contracting a disease or infection.
3. Before occupational use of this respirator, a written respiratory protection program must be implemented meeting all the local government requirements. In the United States employers must comply with OSHA 29 CFR 1910.134 which includes medical evaluation, training, and fit testing.

**Makrite**<sup>®</sup>  
Life Holds Unlimited Possibilities

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## Points For Attention

1. This Report Is Invalid If Without “The Text Report Special Seal” Or The Official Seal Of The Institute.
2. The Reproduction Is Invalid If Without Being Confirmed By “The Text Report Special Seal” Or The Official Seal Of The Institute.
3. This Report Is Invalid If Without Signature Of The Major Tester And The Examiner And The Approver.
4. This Report Is Invalid If In Any Form By Any Means Altered.

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The Institute E-mail:1735889887@qq.com

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# NIOSH certification

DEPARTMENT OF HEALTH & HUMAN SERVICES

NIOSH Reference: TN-20673  
Mfr. Reference: MAK-1603

Centers for Disease Control  
and Prevention (CDC)

National Institute for Occupational  
Safety and Health (NIOSH)  
National Personal Protective  
Technology Laboratory (NPPTL)  
626 Cochran's Mill Road  
Pittsburgh, PA 15236-0070  
Phone: 412-386-4000  
Fax: 412-386-4051

March 24, 2016

Mr. Alexander Freedman  
Makrite Industries, Inc.  
105 Palmer River Road  
Swansea, MA 02777

Dear Mr. Freedman:

The National Institute for Occupational Safety and Health (NIOSH) has reviewed your request accepted February 8, 2016. This request was for approval of the SEKURA-N95 air purifying filtering facepiece respirator for protections against particulates at a N95 filter efficiency level, reference assembly matrix SEKURA-N95AMaR0.xls, dated: 01/19/2016.

This request is granted. Approvals are granted only for documentation written in the English language. It is the manufacturer's responsibility to correctly translate materials desired in languages other than English. Approval number TC-84A-7638 has been assigned. This respirator is approved for protection against particulates at a N95 filter efficiency level.

The final respirator approval label is included as an attachment to this letter. The abbreviated label has been accepted as submitted. The cautions and limitations which apply to this approval are on the approval label. Only those assemblies affected by this request, or where new approval numbers are assigned, apply to this approval action. Production approval labels cannot include information on unapproved configurations.

The approved assembly consists of the parts as listed on the approval label and the assembly matrix. Parts are to be marked with the numbers indicated on the approval label in a legible and permanent manner (marking cannot be removed without evidence of its previous presence).

This certificate of approval is not an endorsement of the respirator by NIOSH, and such endorsement shall not be stated or implied in advertisements or other publicity. However, you may publicize the fact that this respirator has met the requirements of Title 42, *Code of Federal Regulations*, Part 84 (42 CFR 84).

Page 2 – Mr. Alexander Freedman – TN-20673

No changes may be made to any respirators and accompanying documentation without prior written approval of NIOSH. Requests for changes must be submitted to NIOSH and a modification of this approval must be granted before changes are made.

Sincerely yours,

A handwritten signature in cursive script, appearing to read "David Chirdon".

David Chirdon  
Chief, Conformity Verification and  
Standards Development Branch  
National Personal Protective Technology Laboratory

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# CQC certification

## QUALITY MANAGEMENT SYSTEM CERTIFICATE

Certificate No. 00119Q310176R2M/4400

We hereby certify that

**MAKRITE INDUSTRIES INC./DONGGUAN SURE\_ON PLASTIC &  
FIBER PRODUCTS CO.,LTD./DONGGUAN MAKRITE  
INDUSTRIES LTD.**

Unified Social Credit Code: 9144190056665074X0

11F-5, No.79, Sec 1, Xin Tai 5th Road, XiZhi Dist., New Taipei City, Taiwan,R.O.C./ ZhuFo Ling Industrial Area,TangXia Town,DongGuan City,GuangDong,China / No.12, GuangYuan Street,168 Industrial Area, ZhuFo Ling,TangXia Town,DongGuan City,GuanDong,China

by reason of its

**Quality Management System**

has been awarded this certificate for compliance with the standard

**GB/T 19001-2016 / ISO 9001:2015**

The Quality Management System Applies in the following area:

MAKRITE INDUSTRIES INC: Sales of Personal Protective Equipment(Fiber Mask, Filtering Respirator and Half Face Respirator) and Houseware Products(Table Clothes, Placemat, Outdoor Cover, Shower Curtain and Other Related Products).  
DONGGUAN SURE\_ON PLASTIC & FIBER PRODUCTS CO., LTD: Design and Manufacture of the Exported Personal Protective Equipment (Fiber Mask, Filtering Respirator and Half Face Respirator).  
DONGGUAN MAKRITE INDUSTRIES CO.LTD.: Design and Manufacture Within The Scope of Qualification of Personal Protective Equipment(Fiber Mask, Filtering Respirator and Half Face Respirator), Manufacture of Houseware Products(Table Clothes, Placemat, Outdoor Cover, Shower Curtain and Other Related Products).

**Certified since: December 30, 2013 Valid from: December 9, 2019 Valid until: December 29, 2022**

After a surveillance cycle, the certificate is valid only when used together with an Acceptance Notice of Surveillance Audit issued by CQC. Please access [www.cqc.com.cn](http://www.cqc.com.cn) for checking validity of the certificate.

This certificate and its relevant information can query in the website of Certification and Accreditation Administration of the People's Republic of China ( [www.cnca.gov.cn](http://www.cnca.gov.cn) ).

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陆梅  
Signed by: Lu Mei

**CHINA QUALITY CERTIFICATION CENTRE**

Section 9, No.188, Nanshuan(the South Fourth Ring Road) Xilu(West Road), Beijing 100070,China  
<http://www.cqc.com.cn>

A 0049753

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THE INTERNATIONAL CERTIFICATION NETWORK

# CERTIFICATE

*CQC has issued an IQNet recognized certificate that the organization:*

**MAKRITE INDUSTRIES INC./DONGGUAN SURE\_ON PLASTIC & FIBER PRODUCTS CO.,LTD./DONGGUAN MAKRITE INDUSTRIES LTD.**

11F-5, No.79, Sec 1, Xin Tai 5th Road, XiZhi Dist., New Taipei City, Taiwan,R.O.C./ ZhuFo Ling Industrial Area,TangXia Town,DongGuan City,GuangDong,China / No.12, GuangYuan Street,168 Industrial Area, ZhuFo Ling,TangXia Town,DongGuan City,GuanDong,China

*has implemented and maintains a*

**Quality Management System**

*for the following scope:*

MAKRITE INDUSTRIES INC: Sales of Personal Protective Equipment(Fiber Mask, Filtering Respirator and Half Face Respirator) and Houseware Products(Table Clothes, Placemat, Outdoor Cover, Shower Curtain and Other Related Products).

DONGGUAN SURE\_ON PLASTIC & FIBER PRODUCTS CO., LTD: Design and Manufacture of the Exported Personal Protective Equipment (Fiber Mask, Filtering Respirator and Half Face Respirator).

DONGGUAN MAKRITE INDUSTRIES CO.LTD.: Design and Manufacture Within The Scope of Qualification of Personal Protective Equipment(Fiber Mask, Filtering Respirator and Half Face Respirator), Manufacture of Houseware Products(Table Clothes, Placemat, Outdoor Cover, Shower Curtain and Other Related Products).

*which fulfils the requirements of the following standard:*

**GB/T 19001-2016 / ISO 9001:2015**

*Issued on:* December 9, 2019

*Expires on:* December 29, 2022

**Registration Number: CN00119Q310176R2M/4400**

*This attestation is directly linked to the IQNet Partner's original certificate and shall not be used as a stand-alone document*

Alex Stoichitoiu  
President of IQNet

Lu Mei  
President of CQC

**IQNet Partners\*:**

AENOR Spain AFNOR Certification France APCER Portugal CCC Cyprus CISQ Italy  
CQC China CQM China CQS Czech Republic Cro Cert Croatia DQS Holding GmbH Germany FCAV Brazil  
FONDONORMA Venezuela ICONTEC Colombia Inspecta Sertifiointi Oy Finland INTECO Costa Rica  
IRAM Argentina JQA Japan KFQ Korea MIRTEC Greece MSZT Hungary Nemko AS Norway NSAI Ireland  
NYCE-SIGE México PCBC Poland Quality Austria Austria RR Russia SII Israel SIQ Slovenia  
SIRIM QAS International Malaysia SQS Switzerland SRAC Romania TEST St Petersburg Russia TSE Turkey YUQS Serbia  
IQNet is represented in the USA by: AFNOR Certification, CISQ, DQS Holding GmbH and NSAI Inc.

\* The list of IQNet partners is valid at the time of issue of this certificate. Updated information is available under [www.iqnet-certification.com](http://www.iqnet-certification.com)

# LA certification

## Certificate of LA Mark of Personal Protective Equipment

This is to certify that the Personal Protective Equipments (detail of the products name, category, specification, model and mark numbers see appendix of LA Mark Certificate) manufactured by **MAKRITE HUBEI INDUSTRIAL CO., LTD** meet the requirements of LA Mark.

Certificate number: **LA-2019-2010**

Date of issue: **Dec.31,2019**

Date of expiry: **Dec.31,2024**

**LA MARK CENTER OF  
PERSONAL PROTECTIVE EQUIPMENT**