



# Nitrile Powder Free

Examination Gloves

3400 Series

- Provides tensile strength and excellent comfort, while not sacrifising tactile sensitivity.
- Ultra 300 Saver Pack reduces storage space and packaging waste.
- Fingertip textured in blue color.
- Meets ASTM D6319 Standards.
- 300 pcs per box (Size XS-L), 250 pcs per box (Size XL) by weight. 10 boxes per carton.
- · Shelf Life: 3 Years
- Available Size: XS 3405

S 3406

M 3407

L 3408

XL 3409



#### Saver Pack Series Benefits

Maximize Storage Space	Reduce Packaging Waste	Save on Shipping Cost	Reduce Stocking Time
100			<b>888</b>
200			••
300 €	Î		<b>©</b>

Features	Benefits	
Non-Latex	- No latex allergies.	
Powder Free	- Eliminates powder allergies and dry skin irritation.	
Finger-Textured	- Provides excellent grip.	
300pcs Saver Pack	<ul> <li>Triple the quantity of gloves packed in the same box</li> <li>Reduces the usage of packaging materials.</li> </ul>	

Test	Particular	Specification	
Dimension (mm)	X-Small (3405)	70 ± 10	
Palm Width	Small (3406)	80 ± 10	
	Medium (3407)	95 ± 10	
	Large (3408)	111 ± 10	
	X-Large(3409)	120 ± 10	
Dimension (mm)	X-Small (3405)	min 220	
Length	Small (3406)		
	Medium (3407)		
	Large (3408)	min 230	
	X-Large (3409)		
Thickness (mm)	Palm Thickness	min. 0.05	
······	Finger Thickness	min. 0.05	
Tensile Strength (MPa)	Before Aging	min. 14	
	After Aging	min. 14	
Elongation (%)	Before Aging	min. 500	
	After Aging	min. 400	
Powder Residue	Powder Free (mg/glove)	max. 2	
Freedom From Holes	G1/AQL 2.5		
Physical Appearance	Rolled beaded cuff		
	Blue color		
	Fingertip-textured surface		
	Ambidextrous		

Remark: The above data represent the norm of the product which may be changed without notice.

Gloves used for protection against chemotherapy drugs exposure should be selected specifically for the type of drugs being used. User is responsible for determining the safe use of glove in chemotherapy treatment. These gloves are for single and transient use only.

## Tested for use with Chemotherapy Drugs per ASTM D6978

Tested Chemotherapy Drug & Concentration	Minimum Breakthrough Detection Time	
Cisplatin, 1.0 mg/ml	>240 mins	
Cyclophosphamide (Cytoxan), 20.0 mg/ml	>240 mins	
Dacarbazine (DTIC), 10.0 mg/ml	>240 mins	
Doxorubicin Hydrochloride, 2.0 mg/ml	>240 mins	
Etoposide (Toposar), 20.0 mg/ml	>240 mins	
Fluorouracil, 50.0 mg/ml	>240 mins	
Paclitaxel (Taxol), 6.0 mg/ml	>240 mins	
Methotrexate, 25 mg/ml	>240 mins	
Mitomycin C, 0.5 mg/ml	>240 mins	
Vincristine Sulfate, 1.0 mg/ml	>240 mins	
Carmustine (BCNU), 3.3 mg/ml	Not recommended	
Thiotepa, 10.0 mg/ml	Not recommended	

Warning: Do not use with Carmustine and Thiotepa

#### **Packaging Dimension**

Inner	Carton	
245mmL x 128mmW x 88mmH	460mmL x 266mmW x 265mmH	

#### **Intended Application**

Cranberry Evolve®CT Nitrile Powder Free Examination Gloves is a disposable glove product that is only applicable for examination, medical treatments, and handling contaminated medical devices to prevent cross contamination betwen users and patients, and not suitable for surgical use.

### **Handling and Storage**

Store in a dry place at room temperature. Exposed product should be shielded from direct sunlight, intense artificial light, x-ray machines and other sources of ozone.

#### **Hazard Identification**

If allergic reaction occurs, discontinue use immediately and consult a physician. his product is not made with natural rubber latex.