

BLUE NITRILE EXAMINATION GLOVES



FEATURES

- Powder-Free
- Ambidextrous
- Beaded Cuff
- Non-Sterile
- Not Made From Natural Rubber Latex

| SPECIFICATIONS | |
|-------------------------|--------------------|
| MATERIAL: | Nitrile |
| POWDER-FREE: | Yes |
| COLOR: | Blue |
| GLOVE LENGTH: | 245mm / 9.5 inches |
| FINGERTIP SURFACE: | Textured |
| PALM THICKNESS: | 0.08mm / 3.1mil |
| FINGER THICKNESS: | 0.11mm / 4.3mil |
| PROCESS: | ONLINE CHLORINATED |
| FDA 510(K) CLEARED: | Yes |
| INSPECTION LEVEL (FFH): | AQL 1.5 |

RECOMMENDED APPLICATIONS

- Medical Examinations
- Laboratory Research
- Food Handling
- Janitorial Work
- Light Duty Maintenance

| ORD | ERIN | IG I | NF | ORI | ИΑТ | ON |
|-----|------|------|----|-----|-----|----|
| | | | | | | |

| ITEM | SIZE | QUANTITY |
|---------------|------|----------|
| NAOCFT245B-S | S | 100/Box |
| NAOCFT245B-M | M | 100/Box |
| NAOCFT245B-L | L | 100/Box |
| NAOCFT245B-XL | XL | 90/Box |

PACKAGING

XS-L: 100 gloves/box 10 boxes/case

XL: 90 gloves/box 10 boxes/case

CHEMICAL PERMEATION TESTING PER ASTM D6978:

RESULTS SHOWN FOR EACH AND EVERY CHEMICAL LISTED: AVERAGE BREAKTHROUGH DETECTION TIME: >240 MINUTES

CHEMICALS TESTED*

CISPLATIN 1.0 MG/ML (1,000PPM)

CYCLOPHOSPHAMIDE 20.0 MG/ML (20,000PPM)

DACARBAZINE 10.0 MG/ML (10,000PPM)

DOXORUBICIN 2.0 MG/ML (2,000PPM)

ETOPOSIDE 20.0 MG/ML (20,000PPM)
FLUOROURACIL 50.0 MG/ML (50,000PPM)
IFOSFAMIDE 50.0 MG/ML (50,000PPM)
MITOXANTRONE 2.0 MG/ML (2,000PPM)

PACLITAXEL 6.0 MG/ML (6,000PPM)

VINCRISTINE SULFATE1.0 MG/ML (1,000PPM)

FENTANYL CITRATE INJECTION 100mcg/2mL

*Additional data available upon request.

